

High Speed Nx/i System includes:

- Large 70 cm gantry aperture
- Up to 350 mA possible
- +/- 30 degree tilt
- 120 sec. spiral
- Scan time 0.8, 1.0, 1.5, 2.0,3 sec
- Patient Positioning with laser light
- Hi-resolution hilite detector
- 0.7 sec partial scanning

Optional Software features:

Smart Prep

Smart mA

Smart Helical

3D imaging Option software

MIP Option software

Real Time Reformation Software

PAN/ZOOM option software

Navigator



Medinnova Systems Pvt. Ltd.

Formerly Soma Tech Pvt. Ltd
Regd Office: E-114 GIDC – Savli
Manjusar - Savli Road, At. Po.
Manjusar - 391 775, Tal. Savli. Dist., Vadodara, Gujarat - India
Phone: +91-2667-264737
imaging@medinnovasystems.com | www.medinnovasystems.com