

Scan Parameters

- .5, .6, .7, .8, .9, 1, 2, 3, 4 second scans
- 80, 100, 120, 140, kVp
- 10 – 440 mA in 10mA increments
- .63, 1.25, 2.5, 3.75, 5, 7.5 10 mm slice thickness
- 1 second minimum interscan delay
- 4 x 912 No. detection channels
- 25, 50 cm scan fields of view
- 120 sec max scan time
- 160 cm max scan volume
- Four slice acquisition



Display

- 2 No flat screen monitors, one for scan, one for display
- 1280 x 1024 Matrices, Pixels
- Up to 16 slices displayed at once

Image Reconstruction

- Silicon Graphics Octane Computer
- .5 sec recon time per slice
- 512 x 512 recon matrix

Generator/Tube

- 53.2 kW maximum power
- 6.3 MHU anode heat capacity
- 840 KHU / min max anode heat dissipation rate

Table/ Gentry

- 51-99 cm vertical range
- 170 cm longitudinal scannable range
- Up to 50cm SFOV
- Support 400 lb patients accurately
- 70cm gantry apparture
- +/- 30 gantry tilt

Medinnova Systems Pvt. Ltd.

Formerly Soma Tech Pvt. Ltd
Regd Office: E-114 GIDC – Savli
Manjusar - Savli Road, At. Po.
Manjusar - 391 775, Tal. Savli. Dist., Vadodara, Gujarat - India
Phone: +91-2667-264737
imaging@medinnovasystems.com | www.medinnovasystems.com

