

## Features:

HiSpeed Fx/i System includes:

- Large 70 cm gantry aperture
- Up to 350 mA Possible
- Plus or Minus (+/-) 30 degree tilt
- Scan time 1.0, 1.5, 2.0,3 sec
- 120 Sec. Spiral
- 1,2 or 3,5 & 10 mm Slice thickness
- Patient Positioning with Laser light
- Hi-resolution Hilite detector
- 0.7 sec partial scanning.

**Optional Software features:**

- Smart Prep
- Smart mA
- Smart Helical
- 3D Imaging Option Software
- MIP Option Software
- Real Time Reformation Software
- PAN / ZOOM option Software
- Display in Display Option Software
- Navigator



**Medinnova Systems Pvt. Ltd.**

Formerly Soma Tech Pvt. Ltd  
Regd Office: E-114 GIDC – Savli  
Manjusar - Savli Road, At. Po.  
Manjusar - 391 775, Tal. Savli. Dist., Vadodara, Gujarat - India  
Phone: +91-2667-264737  
imaging@medinnovasystems.com | www.medinnovasystems.com