

SIGNA Creator covers an amazing scope of utilizations, is intended for convenience and conveys brilliant picture quality and patient solace while improving work process.



Benefits:

- Its demonstrated, exceptionally homogeneous magnet consistently delivers high quality images across a larger FOV than of up to 27% than over conventional analog signal receivers, improving image quality and clinical confidence.
- SIGNA Creator's "head to toe" portfolio of applications gives you neuro performance with automated motion correction, MSK performance with cartilage mapping, multi-planar body imaging with quantitative capabilities, and vascular performance most other MR systems.
- OpTix Optical RF technology offers high quality analog to digital signal conversion and can provide a gain in SNR with non-contrast and automated time-resolved imaging.
- Express Suite coil configuration accomplishes extraordinary inclusion and signal penetration.
- Automatic coil selection helps enable quick and consistent exams.
- The Slide Bar reduces 30 inputs with a single control, simplifying scan optimization for new technologists and experienced users, meanwhile maintaining full control of all parameters.
- Volumetric imaging acquisitions like Cube replace cumbersome, slice-by-slice, plane-after-plane 2D acquisitions with a single 3D volume scan.
- Movement revision methods like PROPELLER help limit the impacts of movement relics, possibly decreasing the requirement for rescans and the effect of patient development on work process.

Features:

- Express coil technology that reduces handling of bulky coils between exams.
- Slide bar, reduces 30 inputs with a single control
- Two low-height table designs that can be lowered to 49cm -
- Fixed table configuration with Express coil technology that enables easy positioning and high throughput
- Detachable table configuration that improves patient access and productivity
- SV25.0 Applications Suite
- 34% less energy consumption
- Smaller footprint, simpler siting
- Motion Correction with Propeller
- Imaging around metal with IDEAL and MAVRIC SL
- Diffusion Imaging with eDWI

Medinnova Systems Pvt. Ltd.

Formerly Soma Tech Pvt. Ltd
Regd Office: E-114 GIDC – Savli
Manjusar - Savli Road, At. Po.
Manjusar - 391 775, Tal. Savli. Dist., Vadodara, Gujarat - India
Phone: +91-2667-264737
imaging@medinnovasystems.com | www.medinnovasystems.com

