

Features:

- General purpose Surgical Table with pendant style hand control.
- 400 lb lift and articulation capacity with distributed load.
- Radiolucent Top and C-arm Compatibility to meet all imaging requirements.
- Top rotates 180° providing improved access for C-arm and Microscope.



Specifications:

Table Top Height	28" Minimum/45" Maximum
Trendelenburg/Revtrend	30°/30°
Lateral Tilt : Right/Left	30°/30°
Flex/Reflex	40°/90°
Head Section Raise/Lower	15°/90° Reversible
Back Section Raise/Lower	90°/40°
Leg Section Raise/Lower	0°/105°
Table top Length	75"
Table top width	20" with side-rails
Patient weight capacity	400lbs lift & articulated (Distributed Load)
Electrical	120 volt, 60 Hz, 5 amp

Medinnova Systems Pvt. Ltd.

Formerly Soma Tech Pvt. Ltd
Regd Office: E-114 GIDC – Savli
Manjusar - Savli Road, At. Po.
Manjusar - 391 775, Tal. Savli. Dist., Vadodara, Gujarat - India
Phone: +91-2667-264737
sales@medinnovasystems.com | www.medinnovasystems.com



Note: The technical data given in this publication is for general information and are subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.