

Features:

- Exceptional imaging capability
 - Radiolucent abductor bars means unobstructed views of crucial patient anatomy
 - Excellent C-arm interface for lower limbs because the table base is completely out of the way
- Designed for Maximum Procedural Flexibility
 - Easy patient transfer and positioning
 - Easy-to-maneuver, lightweight abductor bars
 - Accommodates supine and lateral positions
 - Simple-to-use numbered and lettered positioning locks
 - Facilitates skeletal or boot traction procedures including hip pinning, Intramedullary (IM) nailing, and tibia/fibula surgery
- Flexible Positioning
- Patients weighing up to 400 lb. (181 kg) can be postured on the OrthoVision table.

