

The Zimmer A.T.S. 1200 Tourniquet System provides accurate pressure monitoring for patient safety and confidence. Microprocessor controls and built in safety features with self-diagnostic tests make monitoring reliable in a single port, dual cuff surgical tourniquet system.



Features:

- Large easy to read display
- Built in safety features
- 45 minute battery backup
- Visual and audio alarms
- Color coded cuffs in a variety of sizes
- Low profile cuffs
- Lightweight unit is easily used on a table or IV pole
- Self-check calibration each time the machine is turned on
- Dual cuffs allow bilateral or Bier Block procedures with critical monitoring
- Cuff lockout – cuff remains inflated even if power is turned off
- Positive locking connectors for leak free inflation

Medinnova Systems Pvt. Ltd.

Formerly Soma Tech Pvt. Ltd
Regd Office: E-114 GIDC – Savli
Manjusar - Savli Road, At. Po.
Manjusar - 391 775, Tal. Savli. Dist., Vadodara, Gujarat - India
Phone: +91-2667-264737
sales@medinnovasystems.com | www.medinnovasystems.com



Note: The technical data given in this publication is for general information and are subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.