

## Features:

- Portable and designed for tabletop or tourniquet stand monitoring
- Automatic, microprocessor-based pneumatic tourniquet system
- Able to independently control and monitor two cuffs in a Bier Block or bilateral procedure



### Medinnova Systems Pvt. Ltd.

Formerly Soma Tech Pvt. Ltd  
Regd Office: E-114 GIDC – Savli  
Manjusar - Savli Road, At. Po.  
Manjusar - 391 775, Tal. Savli. Dist., Vadodara, Gujarat - India  
Phone: +91-2667-264737  
sales@medinnovasystems.com | www.medinnovasystems.com