

Pentax EG 2930/ EG 2901

Optical System	Field Of View	120°
	Direction Of View	Forward
	Depth Of Field	5mm to 100mm
Insertion Tube	IT Diameter	9.8mm
Angulation	Up	210°
	Down	120°
	Right	120°
	Left	120°
Length	Working Length	105cm
	Total Length	137cm
Biopsy Channel	Channel Size	2.8mm



Pentax EG 3400

Optical System	Field Of View	100°
	Direction Of View	Forward
	Depth Of Field	5mm to 100mm
Insertion Tube	IT Diameter	11.3mm
Angulation	Up	210°
	Down	120°
	Right	120°
	Left	120°
Length	Working Length	105cm
	Total Length	137cm
Biopsy Channel	Channel Size	3.8mm



Medinnova Systems Pvt. Ltd.

Formerly Soma Tech Pvt. Ltd
 Regd Office: E-114 GIDC – Savli
 Manjusar - Savli Road, At. Po.
 Manjusar - 391 775, Tal. Savli. Dist., Vadodara, Gujarat - India
 Phone: +91-2667-264737
 sales@medinnovasystems.com | www.medinnovasystems.com