

Features:

- The Karl Storz Electronic CO₂ endoflator is an insufflation device for universal application in laparoscopic examinations and operations.
- The endoflator provides CO₂ gas distention of the abdomen for the diagnostic and/or operative laparoscopy.
- The device is intended to be used with a central gas supply line to achieve pressures from 0-30 mmHg and flow rates of 0-30 L/min.



Medinnova Systems Pvt. Ltd.

Formerly Soma Tech Pvt. Ltd
Regd Office: E-114 GIDC – Savli
Manjusar - Savli Road, At. Po.
Manjusar - 391 775, Tal. Savli. Dist., Vadodara, Gujarat - India
Phone: +91-2667-264737
sales@medinnovasystems.com | www.medinnovasystems.com