

Specification:

DIMENSIONS (HXWXD) CM, (IN) (DISPLAY)	15x 36 x 30 (5.9x 14.2 x 11.8)
WEIGHT, kg (lb) (DISPLAY)	12.5(27.6)
Type (CONFIGURATION)	Electronic
FLOW OPTIONS, L/MIN	
Range (WAVELENGTH, nm)	0
Fixed settings (FLOW OPTIONS, L/MIN)	1,14,20
VIDEO MONITOR CAPABILITY	Yes
BACTERIAL FILTER (FLOW OPTIONS, L/MIN)	Yes
HIGH-PRESSURE SAFEGUARD (EXTERNAL GAS CYLINDER)	Yes
Alert signal	Audible and visual
Activation press, mm hg (EXTERNAL GAS CYLINDER)	30
Time delay, sec (BREATH DETECTOR)	5
Electronic pressure reduction (EXTERNAL GAS CYLINDER)	Yes
Activation press, mm hg (EXTERNAL GAS CYLINDER)	3-5 mm Hg, programmable
Time delay, sec (BREATH DETECTOR)	2-5, programmable
PERITONEAL PRESSURE INDICATOR (EXTERNAL GAS CYLINDER)	Display
INSUFFLATION VOLUME INDICATOR (FLOW OPTIONS, L/MIN)	Display
EXTERNAL GAS CYLINDER	
Gauge/display (EXTERNAL GAS CYLINDER)	Display
Low-pressure warning (EXTERNAL GAS CYLINDER)	Yes
GAS USED (Hinge)	CO2



Medinnova Systems Pvt. Ltd.

Formerly Soma Tech Pvt. Ltd
 Regd Office: E-114 GIDC – Savli
 Manjusar - Savli Road, At. Po.
 Manjusar - 391 775, Tal. Savli. Dist., Vadodara, Gujarat - India
 Phone: +91-2667-264737
 sales@medinnovasystems.com | www.medinnovasystems.com