

Specification:

- Illumination level setting (connections) :** discrete
- Number of steps (microscope) :** 18
- Mode (imaging system) :** automatic
- Power needed, vac (energy) :** 120
- Hz (operating) :** 50/6
- Dimensions (hxwxd) cm, (in) (display) :** 12.5X 29.5 X 37 (5x 11.6 X 14.6)
- Replacement bulb :** 625
- Type (configuration) :** multipurpose
- Attachment (collimators) :** yes
- Automatic light adjustment (video) :** yes
- Equivalent light cables (connections) :** olympus
- High-intensity light (connections) :** yes
- Lamp used (digital capture) :** xenon - for details please see comments.
- Watts (generator) :** 300
- Of lamps :** 1
- Backup lamp(s) :** 1halogen
- Color temp, \hat{a}° k (connections) :** 6000
 - Lamp duration, hr : 300



Medinnova Systems Pvt. Ltd.

Formerly Soma Tech Pvt. Ltd
Regd Office: E-114 GIDC – Savli
Manjusar - Savli Road, At. Po.
Manjusar - 391 775, Tal. Savli. Dist., Vadodara, Gujarat - India
Phone: +91-2667-264737
sales@medinnovasystems.com | www.medinnovasystems.com