

## Features:

- Compatibility
  - Instrument port : Storz mini
  - Light port adaptor : ACMI, Wolf, Storz/Olympus
  - Eyepiece : Ocular according DIN 58105 threaded
- Objective: Tension-free objective composed-objective
- Image Transmission: Optical fiber bundle system for optimal sharpness, colorsync and resolution
- Utilizes new electroluminescent display and infrared touch screen control technology.
- Light-guiding system: High quality light fibers
- Optical glass: Distal and proximal sapphire
- Biocompatibility: All materials that come in contact with the body are biocompatible. All outside metal parts are made of high quality stainless steel



Diameter In mm	Direction of view in degree	Working length in mm
2,7	0	302
2,7	15	302
2,7	30	302
2,7	70	302
2,9	0	302
2,9	15	302
2,9	30	302
2,9	70	302
3,0	0	302
3,0	15	302
3,0	30	302
3,0	70	302
4,0	0	302
4,0	15	302
4,0	30	302
4,0	17	302

- All Autoclavable

**Medinnova Systems Pvt. Ltd.**

Formerly Soma Tech Pvt. Ltd  
 Regd Office: E-114 GIDC – Savli  
 Manjusar - Savli Road, At. Po.  
 Manjusar - 391 775, Tal. Savli, Dist., Vadodara, Gujarat - India  
 Phone: +91-2667-264737  
 sales@medinnovasystems.com | www.medinnovasystems.com

