

Features:

- **Compatibility**
 - Light port adaptor : ACMI, Wolf, Storz/Olympus
 - Eyepiece : Ocular according DIN 58105 threaded
- **Objective: Tension-free objective composed-objective**
- **Image Transmission: CAD Designed rod-lens system for optimal sharpness, colorsync and resolution.**
- **Light-guiding system: High quality light fibers**
- **Optical glass: Distal and proximal sapphire**
- **Biocompatibility: All materials that come in contact with the body are biocompatible. All outside metal parts are made of high quality stainless steel**



Diameter In mm	Direction of view in degree	Working length in mm
2,9	0	302
2,9	30	302
5,0	0	302
5,0	30	302
10	0	316
10	30	316
10	0	349
10	30	349
10	0	431
10	30	431

- All Autoclavable

Medinnova Systems Pvt. Ltd.

Formerly Soma Tech Pvt. Ltd
 Regd Office: E-114 GIDC – Savli
 Manjusar - Savli Road, At. Po.
 Manjusar - 391 775, Tal. Savli. Dist., Vadodara, Gujarat - India
 Phone: +91-2667-264737
 sales@medinnovasystems.com | www.medinnovasystems.com

